



## INSTRUMENT CLEANING AND PRE-DISINFECTION

### Anios'Clean Excel D

#### FOR THE CLEANING AND DISINFECTION OF SURGICAL AND MEDICAL INSTRUMENTS

- Patented concept
- Original combination of active substances in synergy
- Non-corrosive, compatible with a broad range of materials




# Anios'Clean Excel D

For the cleaning and disinfection of surgical and medical instruments


## PRODUCT CHARACTERISTICS

- Clear blue solution
- Pure product pH: min 6.5 - max 7.5
- pH at 0.5% in tap water: 7.4
- pH at 0.5% in soft water: 7.6
- pH at 0.5% in demineralized water: 6.6


## INSTRUCTIONS FOR USE




Concentrated solution for professional use. Use for soaking at a required solution, respecting the minimum contact time from 5 minutes.



Totally immerse the medical device. Cover the basin. Keep the soaking bath covered during the contact time. Brush if visible soiling is present.



Thoroughly rinse treated medical devices with tap water to remove all traces of product and organic soiling.



Replace the solution as soon as it is visibly soiled and at least once per day to maintain the effectiveness of the product. Prime the associated pump before first use.

## DOSING CHART

LITRE SOLUTION	CONCENTRATION OF THE APPLICATION SOLUTION		
	0.5%	1.0%	2.0%
1 litre	5 ml	10 ml	20 ml
2 litre	10 ml	20 ml	40 ml
4 litre	20 ml	40 ml	80 ml
8 litre	40 ml	80 ml	160 ml

## MICROBIOLOGICAL PROPERTIES

CLAIM/CONDITIONS	CONCENTRATION	CONTACT TIME
Bactericidal - According to EN 13727 and EN 14561, dirty conditions	0,5%	5 min
Yeasticidal - According to EN 13624 and EN 14562, dirty conditions	0,5%	15 min
Tuberculocidal - According to EN 14348 and EN 14563, dirty conditions	2,0%	60 min
Effective against enveloped viruses - According to EN 14476 / EN 17111, dirty conditions	0,5% 1,0%	15 min 5 min

## PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Anios'Clean Excel D	200 x 25 ml doses	2416.097
Anios'Clean Excel D	12 x 1 L dosing bottles	2416.095
Anios'Clean Excel D	4 x 5 L cans with 25 ml pump	2416.036
Anios'Clean Excel D	2 x 5 L cans for DG3 system	2416.038



**CE**  
0459

**LABORATOIRES ANIOS**  
1, Rue de l'Espoir  
59260 Lezennes  
France

Distributed by:  
**LABORATOIRES ANIOS**  
1, Rue de l'Espoir  
59260 Lezennes  
France

**ECOLAB**

**LABORATOIRES ANIOS**  
An Ecolab Company

www.ecolab.com

www.anios.com

## APPLICATION AREA

Cleaning and disinfection of medical devices: surgical and medical instruments, heat-sensitive instruments and endoscopy equipment, before sterilization.

## INSTRUCTIONS FOR USE

Concentrated solution for professional use. Use for soaking at a required solution, respecting the minimum contact time from 5 minutes.

Totally immerse the medical device. Cover the basin. Keep the soaking bath covered during the contact time. Brush if visible soiling is present.

Thoroughly rinse treated medical devices with tap water to remove all traces of product and organic soiling.

Replace the solution as soon as it is visibly soiled and at least once per day to maintain the effectiveness of the product. Prime the associated pump before first use.

## INGREDIENTS

Didecylidimethylammonium chloride, chlorhexidine digluconate, non-ionic surfactants, perfume, colouring, excipients.

## STORAGE


Dangerous - respect the precautions for use (drawn up according to the European rules in force regarding the classification and labelling of chemical products). Storage: between +5°C and +25°C.

## CERTIFICATES

Medical device. Full antimicrobial activity for Anios Clean Excel D is included in the scientific dossier, available on request.

## MEDICAL DEVICE USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.