



HIGH LEVEL DISINFECTANT

Opaster'Anios

HIGH LEVEL DISINFECTANT FOR MEDICAL DEVICES

- Ready to use solution
- Broad antimicrobial spectrum
- Control of the bath conformity using a specific test strip
- Stability of soaking bath: 14 days




Opaster'Anios

High level disinfectant for medical devices


PRODUCT CHARACTERISTICS

- Clear blue solution
- Active on Clostridium difficile
- 0.55% orthophthalaldehyde solution buffered at pH 7
- Broad compatibility with all types of materials


INSTRUCTIONS FOR USE



Pre-disinfection step: Clean the medical device with a detergent or a pre-disinfectant product. Rinse thoroughly. For endoscopic equipment: brush channels.




Disinfection step: Pour the solution. Comply with the minimum contact time of 5 minutes. Keep the soaking bath covered during the contact time. For endoscopic equipment: make solution circulate inside channels. Do not add any product into OPASTER'ANIOS.



Rinse thoroughly with sterile or filtered water (0.2 µm) to eliminate all traces of the product.

For endoscopic equipment: rinse the internal and external parts of the medical device.



Dry with a single-use towel. Keep the medical device as aseptically as possible up to use it again. To guarantee the efficacy of the disinfectant, Check the ortho-phthalaldehyd level with an OPASTER'ANIOS testing strip at the start of the session, when restarting after an interruption and at least once a day. Maximum storage of the soaking bath: 14 days (except France - 7 days maximum acc. French regulation).

Check the ortho-phthalaldehyd level with an OPASTER'ANIOS testing strip at the start of the session, when restarting after an interruption and at least once a day.

MICROBIOLOGICAL PROPERTIES

CLAIM/CONDITIONS	CONTACT TIME
Bactericidal - According to EN 13727 / EN 14561, clean conditions	5 min
Yeasticidal - According to EN 13624 / EN 14562, clean conditions	5 min
Fungal - According to EN 13624 / EN 14562, clean conditions	5 min
Tuberculocidal - According to EN 14348/EN 14563, clean conditions	5 min
Mycobactericidal - According to EN 14348 / EN 14563, clean conditions	5 min
Effective against enveloped viruses According to EN 14476/EN 17111, clean conditions	5 min
Virucidal - According to EN 14476 / EN 17111, clean conditions	5 min
Sporicidal against C. Difficile - According to EN 17126, clean conditions	30 min

PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Opaster'Anios	4 x 5 L cans	2504.028
Opaster'Anios	Test strips (box of 50)	100.110



LABORATOIRES ANIOS
1, Rue de l'Espoir
59260 Lezennes
France

Distributed by:
LABORATOIRES ANIOS
1, Rue de l'Espoir
59260 Lezennes
France

ECOLAB

LABORATOIRES ANIOS
An Ecolab Company

www.ecolab.com

www.anios.com

APPLICATION AREA

High level disinfectant for medical devices: surgical and medical equipment, endoscopy equipment and heat-sensitive equipment.

INSTRUCTIONS FOR USE

Ready-to-use solution for professional use. Do not dilute. Carefully clean and rinse the medical device before soaking it in OPASTER'ANIOS.

Pre-disinfection step: Clean the medical device with a detergent or a pre-disinfectant product. Rinse thoroughly. For endoscopic equipment: brush channels.

Disinfection step: Pour the solution. Comply with the minimum contact time of 5 minutes. Keep the soaking bath covered during the contact time. For endoscopic equipment: make solution circulate inside channels. Do not add any product into OPASTER'ANIOS. Rinse thoroughly with sterile or filtered water (0.2 µm) to eliminate all traces of the product.

For endoscopic equipment: rinse the internal and external parts of the medical device. Dry with a single-use towel. Keep the medical device as aseptically as possible up to use it again.

To guarantee the efficacy of the disinfectant, Check the ortho-phthalaldehyd level with an OPASTER'ANIOS testing strip at the start of the session, when restarting after an interruption and at least once a day. Maximum storage of the soaking bath: 14 days (except France - 7 days maximum acc. French regulation).

INGREDIENTS

Ortho-Phthalaldehyde solution (0.55%) in a buffered medium, colouring, perfume.

STORAGE


Dangerous - respect the precautions for use (drawn up according to the European rules in force regarding the classification and labelling of chemical products). Storage: between +5°C and +25°C.

CERTIFICATES

Medical device. Full antimicrobial activity for Opaster'Anios is included in the scientific dossier, available on request.

MEDICAL DEVICE USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.