



READY TO USE HIGH LEVEL DISINFECTANT

Steranios 2%

**FOR MEDICAL DEVICES SURGICAL AND MEDICAL
EQUIPMENT, ENDOSCOPY EQUIPMENT AND
HEAT-SENSITIVE EQUIPMENT**

- Solution ready to use: absence of activator
- Active against bacteria, yeasts, moulds, virus and mycobacteria in 10 minutes
- Possible control of disinfectant solution's conformity with test-strips (every 4 hours)




Steranios 2%

For medical devices surgical and medical equipment, endoscopy equipment and heat-sensitive equipment


PRODUCT CHARACTERISTICS

- Clear green solution
- 2% glutaraldehyde solution buffered at pH 6 (sodium citrate)
- Wide compatibility with all kind of materials


INSTRUCTIONS FOR USE




1 Pre-disinfection step: Clean the medical device with a detergent or a pre-disinfectant product. Rinse thoroughly. For endoscopic equipment: brush channels.



2 Disinfection step: Pour the solution. Keep the soaking bath covered during the contact time. Respect the minimum contact time. For endoscopic equipment: make solution circulate inside channels.



3 Rinse thoroughly with sterile or filtered water (0.2 µm) to eliminate any trace of the product. For endoscopic equipment: rinse the internal and external parts of the medical device.



4 Dry with a single-use towel. Keep the medical device as aseptically as possible up to use it again. Maximum storage of the soaking bath: 30 days.

MICROBIOLOGICAL PROPERTIES

CLAIM/CONDITIONS	CONTACT TIME
Bactericidal - According to EN 13727/EN 14561, clean conditions	5 min
Yeasticidal - According to EN 13624/EN 14562, clean conditions	5 min
Fungicidal - According to EN 13624/EN 14562, clean conditions	10 min
Tuberculocidal - According to EN 14348/EN 1463, clean conditions	10 min
Mycobactericidal - According to EN 14348/EN 14563, clean conditions	10 min
Virucidal - According to EN 14476/EN 17111, clean conditions	10 min
Sporicidal against C. Difficile - According to EN 17126, clean conditions	60 min

PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Steranios 2%	4 x 5 L canister	382.034
Steranios 2%	Testing strips of Steranios 2%	100.987



CE
0459

LABORATOIRES ANIOS
1, Rue de l'Espoir
59260 Lezennes
France

Distributed by:
LABORATOIRES ANIOS
1, Rue de l'Espoir
59260 Lezennes
France

ECOLAB

LABORATOIRES ANIOS
An Ecolab Company

www.ecolab.com

www.anios.com

APPLICATION AREA

High level disinfectant for medical devices surgical and medical equipment, endoscopy equipment and heat-sensitive equipment.

INSTRUCTIONS FOR USE

Ready-to-use product for professional use only. Carefully clean and rinse the medical device before soaking it in STERANIOS 2%.

- 1.** Pre-disinfection step: Clean the medical device with a detergent or a pre-disinfectant product. Rinse thoroughly. For endoscopic equipment: brush channels.
- 2.** Disinfection step: Pour the solution. Keep the soaking bath covered during the contact time. Respect the minimum contact time. For endoscopic equipment: make solution circulate inside channels.
- 3.** Rinse thoroughly with sterile or filtered water (0.2 µm) to eliminate any trace of the product. For endoscopic equipment: rinse the internal and external parts of the medical device.
- 4.** Dry with a single-use towel. Keep the medical device as aseptically as possible up to use it again. Maximum storage of the soaking bath: 30 days.

INGREDIENTS

Steranios 2% is a 2% glutaraldehyde solution (expressed in 100% active compound), buffered at pH 6 in the presence of surface effects catalysor.

STORAGE


Dangerous - respect the precautions for use (drawn up according to the European rules in force regarding the classification and labelling of chemical products).
Storage: from +5°C to +25°C.

CERTIFICATES

Medical device. Full antimicrobial activity for Steranios 2% is included in the scientific dossier, available on request.

MEDICAL DEVICE USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.