

MOBILE MODULE FOR OPERATING THEATRE PLATFORM

Mobile module for collection and pre-disinfection of instrumentation platforms in operating theater.

Instrubac for operating theater platform :

40L soaking bin, full height graduated, equipped with diamond-shaped tank bottom and sump to allow draining.

Dimensions :

Total height with Instrubac : 962 mm
Module : L 850 x W 445 x H 917 mm
Instrubac for operating theater platform : L 730 x W 455 x H 345 mm (external dimensions)

Module technical characteristics :

- Stainless steel
- Equipped with 4 pivoting wheels including 2 with brakes.
- Equipped with fall protection for inferior soaking bin.



Composition	REF
1 Set without draining : trolley + 40L soaking bin + lid	404 660
2 Total draining set : trolley + 40L soaking bin + lid	404 700
3 Quick draining set : trolley + 40L soaking bin + lid	404 880
Accessories	REF
Traceability plate (plate + pliers + timer)	404 661
Curved lid for 40L soaking bin	404 656
Trolley only	404 659
Soaking bin 40 L without draining	404 655
Soaking bin 40 L quick draining	404 655B
Soaking 40 L total draining	404 655A

MOBILE MODULE FOR ENDOSCOPY

Soaking bins not delivered with module



TECHNICAL CHARACTERISTICS

Stainless steel mobile for 15L or 20L INSTRUBACS suitable for autoclaving.

Dimensions

Module for INSTRUBACS 15L :

970 L x 440 W x 760 H mm
1035 L x 505 W x 760 H mm (external dimensions)
Height with Instrubac 15L : 870mm

Module for INSTRUBACS 20L :

1095 L x 350 W x 760 H mm
1200 L x 410 W x 760 H mm (external dimensions)
Height with Instrubac 20L : 920mm

References

Module for Instrubac 15L.....404 077
Module for Instrubac 20L.....404 020

[To order the soaking bins, refer to the leaflet].



TREATMENT OF INSTRUMENTATION

INSTRUBACS RANGE



- INSTRUBACS FOR AUTOCLAVING
- INSTRUBACS NOT SUITABLE FOR AUTOCLAVING
- MOBILE MODULES



Pavé du Moulin - 59260 Lille-Hellemmes - FRANCE

Tél. +33 3 20 67 67 67 - Fax +33 3 20 67 67 68

www.anios.com

INSTRUBACS FOR AUTOCLAVING


SOAKING BINS FOR COLLECTION, PRE-DISINFECTION AND DISINFECTION OF INSTRUMENTATION

SOAKING BIN FOR DRILL



Live vol	External Dim.*	References
0,3 L	Ø 8,5 x H6,8 cm	404 050

1L MICROSURGERY SOAKING BIN



> Lid with hinges

Live vol	External Dim.*	Dim. strainer*	References
0,5 L	L26,5 x W14 x H9	L24,5 x W11 x H7,5 (2,7)	404 010

20L ENDOSCOPY SOAKING BIN MODULE

FOR ENDOSCOPY



Live vol	External Dim.*	References
20L	L57 x W43 x H30	405 173

20L INSTRUBAC



> Lid with hinges

Live vol	External Dim.*	Dim. strainer*	Without tap	With tap
10 L	L51 x W30,5 x H21	L48,5 x W27 x H17	Simple 404 014	404 021
10 L	L51 x W30,5 x H21	L48,5 x W27 x H17	Double 404 019	404 022

Simple - 1 soaking bin with lid and hinges + 1 strainer
Double - 2 soaking bins with lid and hinges + 1 strainer

INSTRUBACS NOT SUITABLE FOR AUTOCLAVING

SOAKING BINS FOR COLLECTION OF SOILED INSTRUMENTATION

2L OR 10L INSTRUBAC



(1) 10 L (2) 2 L

Live vol	External Dim.*	Dim. strainer*	Without crack	With crack
(1) 5 L	L43 x W29 x H12,5	L36,5 x W25,5 x H9	404 002NM	-
(2) 2 L	L29,5 x W20 x H13	L23 x W17 x H9	404 039	404 038

20L INSTRUBAC



> Lid with hinges

Live vol	External Dim.*	Dim. strainer*	Without tap
10 L	L51 x W30,5 x H21	L48,5 x W27 x H17	404 006

Right tap	Central tap	Slip lid	Left tap
404 023	404 023C	404 023D	404 023G

FLAT 20L INSTRUBAC




FOR ENDOSCOPY

Live vol	External Dim.*	Internal Dim. *	Without tap
10 L	L60 x W40 x H13	L55 x W36,5 x H11,5	406 011

Small side tap	Large side tap
404 411	404 411B

40L and 60L INSTRUBAC



> Better mechanical resistance
> No cutting edge
> Made in polypropylene

Live vol	External Dim.*	Internal Dim.*	Without tap
(1) 20 L	L60 x W40 x H24	L51 x W34 x H22	406 012

Small side tap	Large side tap
404 412	404 412A

Live vol	External Dim.*	Internal Dim.*	Without tap
(2) 40 L	L60 x W40 x H35	L51 x W34 x H30,5	406 013

Small side tap
404 413

20L INSTRUBAC

FOR ORTHOPAEDICS



Live vol	External Dim.*	Internal Dim.*	Without tap
10 L	L59 x W31 x H26,5	L54 x W25,5 x H25	404 179

FURNITURES & EQUIPMENT

PROTECTANIOS

Alternative to soaking bin autoclaving.
Sterile slipcover for Instrubacs (15, 20 and 60 L) and endoscopy soaking bin.



Characteristics

Polythene slipcover designed to fit perfectly the shapes of each soaking bin used for rinsing and avoiding autoclaving step.

Packing - References

	Ref. slipcover	Application	Ref. soaking bin	Pack.
1	405 460	Soaking bin 15 L	404 450	60
2	405 461	Soaking bin 20 L	404 006	60
3	405 462	Belt	405 331	60
		Bac module endoscopie	405 173	
4	405 662	Soaking bin 60 L	406 013	50
5	Belt	Protectanios	405 462	unit

DISINFECTION MOBILE MODULE

Mobile module for collection of stained instruments and non-immersible equipment.
Designed for use with Instrubacs 40L or 60L.



TECHNICAL CHARACTERISTICS

Stainless Steel
Dimensions : 735 L x 490 W x 735 H (mm)
External dimensions : 765 L x 550 W x 735 H (mm)
Height with 40L INSTRUBAC: 950mm
Height with 60L INSTRUBAC: 1150mm

Reference


Trolley.....REF. 404 009

*Data expressed in cm

INSTRUBACS NOT SUITABLE FOR AUTOCLAVING

15L INSTRUBAC

FOR ENDOSCOPY



Live vol	External Dim.*	Dim. strainer*	Without tap	With tap
(1) 15 L	L90 x W21 x H16	L85,5 x W18 x H7,5	Simple 404 450	404 024
(2) 15 L	L90 x W21 x H16	L85,5 x W18 x H7,5	Double 404 451	404 025

(1) - 1 soaking bin + 1 strainer + 1 lid + 1 platform
(2) - 2 soaking bins + 1 strainer + 2 lids + 1 platform

4L INSTRUBAC



Live vol	External Dim.*	Dim. strainer*	Simple	Double
4L	L47 x W18 x H12,5	L42 x W14 x H8,5	404 157	404 158

(1) - 1 soaking bin + 1 strainer + 1 lid + 1 platform
(2) - 2 soaking bins + 1 strainer + 2 lids + 1 platform

*Data expressed in cm